

GST-Trap_M for Immunoprecipitation of GST-Fusion Proteins

For the immunoprecipitation of GST-fusion-proteins from cellular extracts.

Only for research applications, not for diagnostic or therapeutic use

1. Introduction For biochemical analyses including mass spectrometry and enzyme activity measurements GST-fusion proteins and their interacting factors can be isolated fast and efficiently in one step by immunoprecipitation using the GST-Trap. GST-Trap_M contains a small GST-binding protein covalently coupled to the surface of magnetic beads. It enables the purification of any protein of interest fused to GST.

2. Content

Reagent	Code	Quantity
GST-Trap_M	stm-20	20 reactions (0.5 ml resin)
GST-Trap_M	stm-100	100 reactions (2.5 ml resin)
GST-Trap_M	stm-200	200 reactions (5 ml resin)
GST-Trap_M	stm-400	400 reactions (10 ml resin)

3. Bead properties Bead size: ~ 0.5 - 1 μm
Storage buffer: 1x PBS, 0.05% sodium azide
Binding capacity: 10 μl GST-Trap_M slurry binds 0.25 – 0.5 μg of GST

4. Stability and Storage Shipped at ambient temperature. Upon receipt store at +4°C.
Stable for 6 month. Do not freeze.

5. Protocol

- For one immunoprecipitation reaction resuspend cell pellet of 1 ml bacteria culture expressing GST in 100 - 200 μl lysis buffer by pipetting.
optional: add 1 mM PMSF and 0.1 mg/ml lysozyme to lysis buffer
- Rotate tube for 1h at 4°C.
- Sonify cell pellet to disrupt bacterial membrane.
- Spin cell lysate at 20.000x g for 5 -10 minutes at 4°C.
- Transfer supernatant to a pre-cooled tube. Adjust volume with dilution buffer to 500 μl . Discard pellet.
For immunoblot analysis dilute 30 μl cell lysate with 10 μl 4x SDS-sample buffer (\rightarrow refer to as input).
- Equilibrate GST-Trap_M beads in dilution buffer. Resuspend magnetic particles by vortexing and transfer calculated volume (20 - 30 μl) in a new reaction cup with 250 μl ice cold dilution buffer. Magnetically separate until supernatant is clear and wash twice with 250 μl of cold dilution buffer.
- Add cell lysate to equilibrated GST-Trap_M beads and incubate the GST-Trap_M beads with the cell lysate under constant mixing for 10 min – 2 h at room temperature or 4°C.
note: during incubation of protein sample with the GST-Trap_M the final concentration of detergents should not exceed 0.2% to avoid unspecific binding to the matrix
- Magnetically separate until supernatant is clear. For western blot analysis dilute 30 μl supernatant with 10 μl 2x SDS-sample buffer (\rightarrow refer to as non-bound). Discard remaining supernatant.

9. Wash beads two times with 500 µl ice cold wash buffer.
optional: increase salt concentration in the second washing step up to 500 mM
10. Resuspend GST-Trap_M beads in 100 µl 2x SDS-Sample buffer or go to step 11.
11. Boil resuspended beads for 10 minutes at 95°C to dissociate the immunocomplexes from the beads. The beads can be collected by centrifugation at 2.500x g for 2 minutes at 4°C and SDS-PAGE is performed with the supernatant (→ refer to as bound).
12. *optional: elute bound proteins by adding 50 µl 0.2 M glycine pH 2.5 (incubation time: 30 sec under constant mixing) followed by centrifugation. Transfer the supernatant to a fresh cup and add 5 µl 1M Tris base (pH 10.4) for neutralization. To increase elution efficiency this step can be repeated.*

Suggested buffer composition

Buffer	Composition
Lysis buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl
Dilution buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl
Wash buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl
Elution buffer	200 mM glycine pH 2.5

Support/ Troubleshooting

Please refer to our FAQ section at www.chromotek.com or contact support@chromotek.com

Related products

GST Toolbox	code
GST-Trap_M Kit	stmk-20
Blocked magnetic beads	bmp-20
GST antibody	6g9

Limited Use Label License

The purchase of this product conveys to the buyer the limited, non-transferable right to use the purchased amount of the product and components of the product to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly or by implication. For information on obtaining additional rights, please contact licensing@chromotek.com.