

# Spot-Label for Immunofluorescence of Spot-Tag<sup>®</sup> Fusion Proteins

*Only for research applications, not for diagnostic or therapeutic use*

**1. Introduction** Small peptide tags are useful for the labelling and detection of proteins using immunostaining, immunoblotting or immunoprecipitation techniques. The ChromoTek Spot-Tag is a short 12 amino acid affinity tag PDRVRAVSHWSS, which can be cloned either N- or C-terminally to a protein of interest. This tag can be efficiently immunostained with the novel Spot-Label affinity reagent. The Spot-Label consists of a small recombinant bivalent alpaca single-domain antibody fragment covalently coupled to a fluorescent dye. Due to its small size, immunostaining of the Spot-Tag with the Spot-Label minimizes the "linkage error" for super-resolution microscopy applications (e.g. STED and dSTORM). In addition, the Spot-Label has a superior tissue penetration rate, better access to the Spot epitope, and higher labelling density.

## 2. Content

Reagent	Quantity	Code
Spot-Label ATTO594 for Immunofluorescence (IF), bivalent	50 µl	eba594-50
Spot-Label ATTO594 for Immunofluorescence (IF), bivalent	10 µl	eba594-10

Concentration: 1 g/L. Storage buffer: 1x PBS, 0.09% sodium azide.

## 3. Optical Properties

**ATTO 594:** Excitation range 580 - 615 nm ( $\lambda_{abs}$ = 601 nm)  
Emission range 620 - 660 nm ( $\lambda_{fl}$ = 627 nm)

For further information please refer to <http://www.atto-tec.com>

## 4. Stability and Storage

Shipped at ambient temperature. Upon receipt store at +4°C / 40°F.  
Stable for 6 months. Do not freeze. Protect from light.

## 5. IF Protocol

- Fixation:** Fix cells seeded on coverslips in 3.7% formaldehyde in PBS for 10 min at room temperature.  
*Note: Always prepare a fresh formaldehyde dilution.*
- Wash samples three times with PBS (Phosphate Buffered Saline). Do not store fixed cells.
- Permeabilization:** Add PBS containing 0.5% Triton X-100 to samples and incubate for 5 min at room temperature.  
*Note: Alternatively, use ice-cold 100% methanol for permeabilization.*
- Wash samples twice with PBS.
- Blocking:** Add 4% BSA in PBS to samples and incubate for 20 min at room temperature.  
*Note: If necessary, use additional blocking reagents (e.g. 10% normal serum in PBS or Image-iT<sup>™</sup> FX Signal Enhancer from ThermoFischer Scientific) and extend the blocking time up to 60 min.*
- Spot-Label incubation:** Dilute Spot-Label 1:2000 in blocking buffer and incubate for overnight at +4°C.  
*Note: For multiplexing protocols, you can combine Spot-Label with another primary or secondary antibody.*
- Wash samples three times for 5-10 min in PBS.
- If required, counterstain with DNA fluorescent dyes, e.g. DAPI in PBS. Proceed with imaging directly or mount samples, if necessary.
- Mounting:** Rinse sample briefly in water to prevent salt crystal formation. Mount in ProLong<sup>™</sup> Diamond Antifade Mountant from ThermoFischer Scientific or other mounting media with anti-fading agents.

### Suggested buffer composition

Buffer	Composition
Fixation buffer	3.7% formaldehyde in PBS
Permeabilization buffer	PBS; 0.5% Triton X-100
Wash buffer	PBS
Blocking buffer	4% BSA (w/v); PBS

### Support/ Troubleshooting

Please refer to our FAQ section at [www.chromotek.com](http://www.chromotek.com) or contact [support@chromotek.com](mailto:support@chromotek.com)

### Related Products

Spot-Tag Toolbox	Code
Spot-Trap® Agarose	eta-20
Spot-Trap® Magnetic Agarose	etma-20
Blocked agarose beads	bab-20
Spot-Label, uncoupled, for WB & IP, monovalent	etx-10; etx-250
Spot-Tag peptide	ep-1
Spin columns	sct-10; sct-20; sct-50

### Trademarks and Copyrights

Spot-Label, Spot-Tag and Spot-Trap are registered trademarks of ChromoTek GmbH. *Image-iT™* and *ProLong™* are trademarks of Thermo Fischer Scientific Inc.