

Myc-Trap[®]_A Kit for immunoprecipitation of Myc-tagged proteins from mammalian cell extract

Only for research applications, not for diagnostic or therapeutic use.

Introduction Epitope tags are useful for the labelling and detection of proteins using immunoblotting, immunoprecipitation, and immunostaining techniques. Because of their small size, they are unlikely to affect the tagged protein's biochemical properties. The Myc epitope tag is widely used to detect expression of recombinant proteins in bacteria, yeast, insect and mammalian cell systems. For biochemical analysis including mass spectrometry and enzyme activity measurements these Myc-tagged proteins and their interacting factors can be isolated fast and efficiently by immunoprecipitation using the Myc-Trap[®]. Myc-Trap[®] utilizes small recombinant antibody fragments covalently coupled to the surface of agarose beads.

Specificity Myc-Trap recognizes Myc-Tag sequence EQKLISEEDL at the N-terminus, C-terminus, or internal site of the fusion protein.

Reagent	Code	Quantity
Myc-Trap [®] _A kit	ytmak-20	20 reactions (0.5 ml resin)
Lysis buffer (CoIP)		30 ml
RIPA buffer		30 ml
5x Wash / Dilution buffer		20 ml

Bead properties Bead size: ~ 90 µm (cross-linked 4% agarose beads)
Storage buffer: 20% EtOH

Stability and Storage Shipped at ambient temperature. Upon receipt store at +4°C.
Stable for 1 year. Do not freeze.

Required solutions **Buffer composition (as provided in the kit)**

Buffer	Composition
Lysis buffer (CoIP)	10 mM Tris/Cl pH 7.5; 150 mM NaCl; 0.5 mM EDTA; 0.5% NP-40, 0.02% Thimerosal
RIPA buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl; 0,5 mM EDTA; 0,1% SDS; 1% Triton X-100; 1% Deoxycholate, 0.02% Thimerosal
5x Dilution/Wash buffer	50 mM Tris/Cl pH 7.5; 750 mM NaCl; 2.5 mM EDTA, 0.02% Thimerosal

Buffer composition (not provided in the kit)

Buffer	Composition
Myc peptide	reconstitute Myc peptide (code yp-1) in water to a final concentration of 2 mg/ ml
2 x SDS-sample buffer	120 mM Tris/Cl pH 6.8; 20% glycerol; 4% SDS, 0.04% bromophenol blue; 10% β-mercaptoethanol

Myc Toolbox	Code
Myc -Trap [®] _A	yta-20
Blocked agarose beads	bab-20
Myc-tag antibody [9E1]	9e1
Myc peptide	yp-1
Spin columns	sct-10; sct-20; sct-50

Harvest cells

For one immunoprecipitation reaction the use of $\sim 10^6 - 10^7$ mammalian cells (approx. one 10-cm dish) expressing a protein of interest is recommended. To harvest adherent cells, aspirate growth medium, add 1 ml ice-cold PBS to cells and scrape cells from dish. Transfer cells to a pre-cooled tube, spin at 500 g for 3 min at +4°C and discard supernatant. Wash cell pellet twice with ice-cold PBS, gently resuspending the cells. After washing:

Lyse cells

1. Resuspend cell pellet in 200 μ l ice-cold lysis buffer by pipetting or using a syringe.
note: Supplement lysis buffer with protease inhibitors and 1 mM PMSF (not included).
optional for nuclear/chromatin proteins: Use RIPA buffer supplemented with 1 mg/ml DNase, 2.5 mM MgCl₂, protease inhibitors and 1 mM PMSF (not included).
2. Place the tube on ice for 30 min with extensively pipetting every 10 min.
3. Centrifuge cell lysate at 20.000x g for 10 min at +4°C. Transfer supernatant to a pre-cooled tube. Add 300 μ l dilution buffer to supernatant. Discard pellet.
note: At this point cell lysate may be put at -80°C for long-term storage.
optional: Add 1 mM PMSF and protease inhibitors (not included) to dilution buffer.

Equilibrate beads

4. Vortex Myc-Trap[®]_A beads and pipette 25 μ l bead slurry into 500 μ l ice-cold dilution buffer. Centrifuge at 2.500x g for 2 min at +4°C. Discard supernatant and repeat wash twice.

Bind proteins

5. Add diluted supernatant (step 3) to equilibrated Myc-Trap[®]_A beads (step 4). If required, save 50 μ l of diluted lysate for immunoblot analysis. Tumble end-over-end for 1 hour at 4°C.
6. Centrifuge at 2.500x g for 2 min at +4°C. If required, save 50 μ l supernatant for immunoblot analysis. Discard remaining supernatant.

Wash beads

7. Resuspend Myc-Trap[®]_A beads in 500 μ l ice-cold dilution buffer. Centrifuge at 2.500x g for 2 min at +4°C. Discard supernatant and repeat wash twice.
optional: Increase salt concentration in the second washing step up to 500 mM.

Elute proteins

8. Resuspend Myc-Trap[®]_A beads in 100 μ l 2x SDS-sample buffer.
9. Boil resuspended Myc-Trap[®]_A beads for 10 min at 95°C to dissociate immunocomplexes from Myc-Trap[®]_A beads. Myc-Trap[®]_A beads can be collected by centrifugation at 2.500x g for 2 min at 4°C and SDS-PAGE is performed with the supernatant.
10. *Optional instead of steps 8 and 9: Reconstitute Myc peptide (code yp-1) in water to a final concentration of 2 mg/ml. Add $\sim 25 \mu$ g Myc peptide in 50-100 μ l in dilution buffer to Myc-Trap[®]_A beads. Mix for 15 min at RT. Centrifuge at 2500x g for 2 min and carefully transfer supernatant to a new tube. To increase elution efficiency this step can be repeated.*
Note: Use our spin column protocol for easy elution

Support

Please refer to our FAQ section at www.chromotek.com or contact support@chromotek.com

Limited Use Label License

The purchase of this product conveys to the buyer the limited, non-transferable right to use the purchased amount of the product and components of the product to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly or by implication. For information on obtaining additional rights, please contact licensing@chromotek.com.

Optional: Spin column protocol for immunoprecipitation and elution of proteins from Myc-Trap[®] agarose beads

Note: spin columns (product code sct-10) are not included

Harvest cells	For one immunoprecipitation reaction the use of $\sim 10^6$ - 10^7 mammalian cells (approx. one 10-cm dish) expressing a protein of interest is recommended. To harvest adherent cells, aspirate growth medium, add 1 ml ice-cold PBS to cells and scrape cells from dish. Transfer cells to a pre-cooled tube, spin at 500 g for 3 min at +4°C and discard supernatant. Wash cell pellet twice with ice-cold PBS, gently resuspending the cells. After washing:
Lyse cells	<ol style="list-style-type: none">1. Resuspend cell pellet in 200 μl ice-cold lysis buffer by pipetting or using a syringe. <i>note: Supplement lysis buffer with protease inhibitors and 1 mM PMSF (not included).</i> <i>optional for nuclear/chromatin proteins: Use RIPA buffer supplemented with 1 mg/ml DNase, 2.5 mM MgCl₂, protease inhibitors and 1 mM PMSF (not included).</i>2. Place the tube on ice for 30 min with extensively pipetting every 10 min.3. Centrifuge cell lysate at 20.000x g for 10 min at +4°C. Transfer supernatant to a pre-cooled tube. Add 300 μl dilution buffer to supernatant. Discard pellet. <i>note: At this point cell lysate may be put at -80°C for long-term storage.</i> <i>optional: Add 1 mM PMSF and protease inhibitors (not included) to dilution buffer.</i>
Equilibrate beads	<ol style="list-style-type: none">4. Remove the upper cap of a new spin column and snap off the plug from the bottom of the spin column (Keep cap and plug!). Place the spin column in a 2 ml tube. Vortex Myc-Trap[®]_A beads and pipette 25 μl bead slurry into the spin column. Immediately add 500 μl ice-cold dilution buffer. Centrifuge at 100x g for 5-10 sec. Discard flow-through and repeat wash twice. Close column with the bottom plug.
Bind proteins	<ol style="list-style-type: none">5. Add diluted supernatant (step 3) to equilibrated Myc-Trap[®]_A beads (step 4). If required, save 50 μl of diluted lysate for immunoblot analysis. Secure the top of the spin column. Tumble end-over-end for 1 hour at 4°C.6. Remove the bottom cap from the spin column and place it in a new 2 ml tube. Centrifuge at 100x g for 5-10 sec. If required, save 50 μl flow-through for immunoblot analysis. Discard remaining flow-through.
Wash beads	<ol style="list-style-type: none">7. Resuspend Myc-Trap[®]_A beads in 500 μl ice-cold dilution buffer. Place spin column in a 2 ml tube and centrifuge at 100x g for 5-10 sec. Discard flow-through and repeat wash twice. Close column with the bottom plug. <i>optional: Increase salt concentration in the second washing step up to 500 mM.</i>
Elute proteins	<ol style="list-style-type: none">8. Reconstitute Myc peptide (code yp-1) in water to a final concentration of 2 mg/ml. Add $\sim 25\mu$g Myc peptide in 50-100μl dilution buffer to Myc-Trap[®]_A beads. Mix for 15 min at RT. Remove bottom plug of the spin column and place it in a new tube. Centrifuge at 1000x g for 30-60 sec. To increase elution efficiency this step can be repeated.

Support Please refer to our FAQ section at www.chromotek.com or contact support@chromotek.com

Limited Use Label License

The purchase of this product conveys to the buyer the limited, non-transferable right to use the purchased amount of the product and components of the product to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly or by implication. For information on obtaining additional rights, please contact licensing@chromotek.com.