

p53 C-term-Trap_A for Immunoprecipitation of p53 from mammalian cell extract

Only for research applications, not for diagnostic or therapeutic use.

Introduction The tumor suppressor protein p53 plays a major role in cell cycle control, apoptosis and senescence. p53-Trap is excellent for fast and efficient one-step isolation of p53 and its interacting factors from cellular extract. Isolated p53 protein may be used further for immunoblot analysis or mass spectrometry. p53-Trap utilizes small recombinant alpaca antibody fragments covalently coupled to the surface of agarose beads.

Specificity Species-Reactivity: tested on human, mouse
Epitope: 302-393

Content	Reagent	Code	Quantity
	p53 C-term-Trap_A	pt2a-20	20 reactions (0.5 ml slurry)
	p53 C-term-Trap_A	pt2a -100	100 reactions (2.5 ml slurry)
	p53 C-term-Trap_A	pt2a -200	200 reactions (5 ml slurry)
	p53 C-term-Trap_A	pt2a -400	400 reactions (10 ml slurry)

Bead properties Bead size: ~ 90 µm (cross-linked 4% agarose beads)
Storage buffer: 20% EtOH

Stability and Storage Shipped at ambient temperature. Upon receipt store at +4°C.
Stable for 1 year. Do not freeze.

Required solutions **Suggested buffer composition**

Buffer	Composition
Lysis buffer (CoIP)	10 mM Tris/Cl pH 7.5; 150 mM NaCl; 0.5 mM EDTA; 0.5% NP-40
RIPA buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl; 0.5 mM EDTA; 0,1% SDS; 1% Triton X-100; 1% Deoxycholate
Dilution/Wash buffer	10 mM Tris/Cl pH 7.5; 150 mM NaCl; 0.5 mM EDTA
Glycine-elution buffer	200 mM glycine pH 2.5
2 x SDS-sample buffer	120 mM Tris/Cl pH 6.8; 20% glycerol; 4% SDS, 0.04% bromophenol blue; 10% β-mercaptoethanol

Related products	p53 Toolbox	Code
	p53-Binding protein C-term	pt2-250
	p53 C-term-Trap_A Kit	pt2ak-20
	p53-Binding protein N-term	pt-250
	p53-Trap_A N-term	pta-20
	Blocked agarose beads	bab-20
	F2H®-Kit p53-Mdm2/4	f2hk-p53/mdm2 f2hk-p53/mdm4

Support Please refer to our FAQ section at www.chromotek.com or contact support@chromotek.com

Protocol for Immunoprecipitation of p53-Fusion Proteins using p53-Trap_A

Harvest cells

For one immunoprecipitation reaction the use of $\sim 10^6$ - 10^7 mammalian cells (approx. one 10-cm dish) is recommended. To harvest adherent cells, aspirate growth medium, add 1 ml ice-cold PBS to cells and scrape cells from dish. Transfer cells to a pre-cooled tube, spin at 500 g for 3 min at +4°C and discard supernatant. Wash cell pellet twice with ice-cold PBS, gently resuspending the cells. After washing:

Lyse cells

1. Resuspend cell pellet in 200 μ l ice-cold lysis buffer by pipetting or using a syringe.
note: Supplement lysis buffer with protease inhibitors and 1 mM PMSF (not included).
optional for nuclear/chromatin proteins: Use RIPA buffer supplemented with 1 mg/ml DNase, 2.5 mM MgCl₂, protease inhibitors and 1 mM PMSF (not included).
2. Place the tube on ice for 30 min with extensively pipetting every 10 min.
3. Centrifuge cell lysate at 20.000x g for 10 min at +4°C. Transfer lysate to a pre-cooled tube. Add 300 μ l dilution buffer to lysate. Discard pellet.
note: At this point cell lysate may be put at -80°C for long-term storage.
optional: Add 1 mM PMSF and protease inhibitors (not included) to dilution buffer.

Equilibrate beads

4. Vortex p53_A beads and pipette 25 μ l bead slurry into 500 μ l ice-cold dilution buffer. Centrifuge at 2.500x g for 2 min at +4°C. Discard supernatant and repeat wash twice.

Bind proteins

5. Add diluted lysate (step 3) to equilibrated p53_A beads (step 4). If required, save 50 μ l of diluted lysate for immunoblot analysis. Tumble end-over-end for 1 hour at 4°C.
6. Centrifuge at 2.500x g for 2 min at +4°C. If required, save 50 μ l supernatant for immunoblot analysis. Discard remaining supernatant.

Wash beads

7. Resuspend p53_A beads in 500 μ l ice-cold dilution buffer. Centrifuge at 2.500x g for 2 min at +4°C. Discard supernatant and repeat wash twice.
optional: Increase salt concentration in the second washing step up to 500 mM.

Elute proteins

8. Resuspend p53_A beads in 100 μ l 2x SDS-sample buffer.
9. Boil resuspended p53_A beads for 10 min at 95°C to dissociate immunocomplexes from p53_A beads. p53_A beads can be collected by centrifugation at 2.500x g for 2 min at 4°C and SDS-PAGE is performed with the supernatant.

Limited Use Label License

The purchase of this product conveys to the buyer the limited, non-transferable right to use the purchased amount of the product and components of the product to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly or by implication. For information on obtaining additional rights, please contact licensing@chromotek.com.