

**Description:** PNGase F PRIME™ is a recombinant glycosidase, cloned from *Flavobacterium meningosepticum*, which catalyzes the cleavage of N-linked oligosaccharides from proteins.

**Biological Source:** *E. coli*.

**Concentration:** 10<sup>6</sup> Units/mL [2.0 mg/mL]

**Molecular Weight:** PNGase F PRIME™ has a molecular weight of approximately 36kDa.

**Physical Form:** PNGase F PRIME™ is supplied as a liquid in 1X PBS (Phosphate Buffered Saline; 137 mM NaCl, 10mM Phosphate, 2.7 mM KCl, pH7.4) at a concentration of >10<sup>6</sup> Units/mL.

**Storage Conditions:** Store at temperatures ranging from +2° to -20°C. Avoid multiple freeze-thaw cycles.

**Unit Definition Assay:** Denatured RNase B (10µg) is incubated with PNGase F PRIME™ for 30 minutes at 37°C and then analyzed by SDS-PAGE. Fully glycosylated RNase B migrates at approximately 17kDa. Deglycosylation is assessed by the presence of deglycosylated RNase B with an apparent molecular weight of 13.7 kDa following staining via Coomassie Brilliant Blue™.

**High-End Criteria:** PNGase F PRIME™ is also designed for use in high-end applications. We therefore include the following rigorous quality release criteria using HPLC/UPLC and Mass Spectrometry Imaging of tissue samples:

- Denatured human IgG (10µg) is incubated with PNGase F PRIME™ for one hour before glycan is labeled with the Waters RapiFluor-MS dye and analyzed by normal phase hydrophilic interaction chromatography (HILIC).
- PNGase F PRIME™ is used for imaging of glycans from tissue sections as described in [Powers et al., PLoS One, 2014, 9(9); p. e106255.] using both a Bruker Daltonics SolariX™ 7T Hybrid FTMS System and a Bruker Daltonics rapifleXTM MALDI Tissue typer.

**Purity:** ≥95% as determined by SDS-PAGE analysis and staining with Coomassie Brilliant Blue™.

### Protein Deglycosylation Using Recombinant PNGase F PRIME™

**Note:** The following protocols are intended as a general guide for protein deglycosylation and may require modification for different glycoprotein substrates.

- Recombinant PNGase F PRIME™ has been optimized for use under non-denaturing conditions, but like many enzyme reactions, activity is substrate

dependent and specific conditions should be determined empirically for each target.

- PNGase F PRIME™ is also fully active when used in denaturing reactions.

### 1. General Protocol for the deglycosylation of proteins under Denaturing Conditions:

**Note:** Deglycosylation of glycoproteins may be visualized by gel-shift on SDS-PAGE, with the deglycosylated product running faster than the glycosylated substrate.

#### ***Materials to Be Supplied By the User:***

- 5% SDS
  - 1M DTT
  - 1X Phosphate Buffered Saline (PBS) (pH 7.4)
  - 10% NP-40
  - Wet Ice
- a) Add up to 50µg of the target glycoprotein in 1X PBS to a final volume of 11µL.
  - b) Add 1µL 5% SDS.
  - c) Add 1µL of 1M DTT.
  - d) Denature sample by heating at 95°C for 10 minutes.
  - e) Cool sample by placing sample on Ice.

**Note:** Other buffers can be used if they are within the acceptable pH range for PNGase F PRIME™, pH 6–10.

- g) Add 2µL of 10% NP-40.
- h) Add 1µL of recombinant PNGase F PRIME™.
- i) Incubate at 37°C for 30 minutes.

### 2. General Protocol for the deglycosylation of proteins under Non-Denaturing Conditions:

**Note:** Deglycosylation under non-denaturing conditions may require increasing both the amount of PNGase F PRIME™ used and the incubation time.

#### ***Materials to Be Supplied By the User:***

- 1X Phosphate Buffered Saline (PBS) (pH 7.4)
- a) Add up to 20µg of glycoprotein in 1X PBS to a final volume of 18µL.
  - b) Add 2µL of recombinant PNGase F PRIME™.
  - c) Incubate at 37°C for 0.5–24 hours.