


● [Brain Slice] Electroporation-mediated gene transfer system applied to cultured CNS neurons




Schematic representation of an electroporation set-up.

A fragment of the mouse embryonic hippocampus was placed on a Millipore membrane filter and 5 $\mu$ l EP buffer containing 1mg/ml of plasmid DNA was applied onto the tissue.

A tungsten needle was attached to the surface of a droplet.

After application of square pulses the tissue fragment was returned to a petri dish containing ice-cold HBSS solution.



Electroporation-mediated expression of fluorescent proteins in hippocampal neurons.

(a-c) Organ culture of hippocampal tissue fragments three days after electroporation with CAG-eGFP (a), T $\alpha$ 1X4-eGFP (b), and  $\beta$ -actin-eGFP (c) expression constructs.

(d, e) A mature hippocampal neuron maintained 14 days in dissociated culture after electroporation of a  $\beta$ -actin-eGFP expression construct. Higher magnification view of the region marked by a rectangle in (d) reveals dendritic spines on the surface of dendritic shafts (arrows in e).

(f, g) A hippocampal neuron 7 days after electroporation of 1:1 mixture of T $\alpha$ 1X4-eGFP and T $\alpha$ 1X4-mRFP1. Both eGFP fluorescence (f) and mRFP1 fluorescence (g) can be observed in a single cell.

(h) Relative fluorescence intensity of hippocampal tissue fragments after electroporation of eGFP-expression plasmids with four different promoter sequences. The tissue fragments were maintained in culture for 4 days, fixed and observed under a confocal microscope. Fluorescence intensities per unit area of the tissue fragments were determined.

(i) 2Relative fluorescence intensity of hippocampal tissue fragments isolated at two different developmental stages and electroporated with  $\beta$ -actin-eGFP. Tissue fragments were maintained for 4 days in culture and subsequently fixed. Fluorescence intensities were measured using a confocal microscope.

Bars = 50  $\mu$ m (a-d, f, g); 10  $\mu$ m (e).

Shigeo Okabe, Department of Anatomy and Cell Biology, School of Medicine, Tokyo Medical and Dental University;

\*Neuroreport, Volume 15, Issue 6, Pages 971-975, April 29, 2004

● [Brain] Electroporation-mediated gene transfer in the adult rat brain




Figure A: EGFP expression in the medial preoptic nuclei of a female rat examined 4 days after bilateral electroporation at 10 weeks of age. (An asterisk indicates the trace of the positioning of the electrode)




Figure B: EGFP-positive cells (high magnification of Fig. A using a 60x objective lens). EGFP fluorescent signals are observed in the perikarya.



Figure C: Estrogen receptor immunoreactivity in the medial preoptic nuclei and the periventricular nuclei of an adult female rat. 3V: third ventricle




Figure D: Estrogen receptor  $\alpha$ -positive cells (high magnification of Fig. C using a 60x objective lens). Estrogen receptor  $\alpha$  immunoreactivity is prominent in the nuclei.

Tetsuo Shirakawa, Center for Advanced Oral Medicine, Hokkaido University Hospital