



Midori Green Advance Tablets

Cat. No.	Pack Size
AG09	75 pcs

Lot no: MG19092017

Exp. Date: 10.2019

Storage: Store at RT, protected from light, shipping at room temperature.

Applications:

- Ideal for routine DNA and RNA gel electrophoresis and blotting assays
- Used dye is non-carcinogenic alternative to Ethidium bromide.

Description: Midori Green Advance Tablets contain everything necessary for an easy preparation of an agarose gel in desired gel percentage.

Midori Green Advance Tablets are packed in a convenient blister pack.

Midori Green Advance Tablets contain:

- Agarose
- TBE powder
- Midori Green Advance DNA stain

Only RT or cooler pure water is needed to add for preparation of an agarose gel.

Do **NOT** use hot water for dissolving the tablet

Do **NOT** add any buffer

This composition is optimized to yield high resolution of sharp DNA bands with high sensitivity and low background.

Midori Green Advance DNA Stain emits green fluorescence when bound to DNA or RNA. It has two secondary fluorescence excitation peaks (~270 nm; ~290 nm) and one strong excitation peak centered around 490 nm. The fluorescence emission is centered at ~530 nm. Thus, Midori Green Advance DNA Stain is compatible with a wide variety of gel reading instruments.

Nippon Genetics Europe GmbH
Binsfelder Straße 77; 52351 Düren; Germany

fax: ++4924215549611, phone: ++492421554960

e-mail: info@nippongenetics.eu

www.nippongenetics.eu/www.nippongenetics.de

Safety: Caution when using hot, viscous solutions! Use suitable safety gear and open bottle gently to avoid accidents.

Midori Green Advance DNA Stain is non-carcinogenic and according to the Ames test it causes significantly fewer mutations than Ethidium bromide. All working solutions can be disposed in regular trash.

Protocol:

- Use the bottle or flask that is at least 3 times of the volume of the solution being prepared.
- Add an appropriate number of agarose tablets in the **water** and do **NOT** add any buffer! See the table below to achieve needed gel percentage.

Gel %	1 tablet	2 tablets
1%	50 ml water	100 ml water
1,5%	33 ml water	67 ml water

- Soak the tablet in the pure **water** for 3-5 minutes (or until it is dissolved) before heating.
- For tablet dissolving use water which is at room temperature, **DO NOT** use hot water.
- Heat the solution until it is clear and visually all the particles are dissolved.
- Cool the gel to 60-70°C and cast the gel, into the gel tray.
- The thickness of gel should be **<0.5cm**.
- Run the gel in TBE running buffer.
- Detect the bands under Blue light or UV illuminator.

Notes:

- Midori Green Advance DNA Stain is non-carcinogenic but may irritate skin and eyes. Please wear gloves while handling.