

Genetics

NEWS

Genetics NIPPON Genetics EUROPE

www.nippongenetics.eu

Products for PCR

efficiency

polymerase

convenience

endpoint PCR

incl. loading dye

direct PCR *incl. dNTPs*

proof-reading

engineered enzyme

archeal type B polymerase

ReadyMix

high GC content

cDNA

FastGene[®] Optima

tissue best results **HotStart**

low copy

long range

SNP **genotyping**

gene knock out

routine PCR

multiplex PCR

mouse tail **very pure**

optimized blend

high-throughput

complex templates

N-T-I-B-O-D-Y

FastGene[®] Optima
(#LS28, 250 units incl. dNTPs)

FastGene[®] Optima HotStart ReadyMix
(#LS29, 500 x 25 µl reactions)

Distributed by:

Bulldog
Bio



Tel.: 603-570-4248
Fax: 603-766-0524
E-Mail: info@bulldog-bio.com
Web: www.bulldog-bio.com

Genotyping of knock-out mice


User's comment

„PCR is essential for us to detect gene-knock out. The FastGene® Optima^{HotStart ReadyMix} suppressed non-specific amplification and produced easy to analyze results. Inexperienced freshman students of the medical school could successfully perform genotyping without making any mistakes.“

The following data were provided by the courtesy of Dr. Mamoru Aoto
Department of Circulatory Physiology, Graduate School of Medicine, Ehime University, Japan.

Experimental conditions

Extraction of DNA from a 2 mm mouse tail tissue:



PCR Reagents:

1. FastGene® Optima HotStart ReadyMix with loading dye (#LS29)
2. Competitor T's enzyme ready mix with loading dye

Template:

1 µl of DNA isolated from mouse tail tissue

PCR product size:

~163 bp

PCR Setup:

FastGene® Optima ^{HotStart ReadyMix}	
Primer set (10 µM)	2.5 µl
Template: mouse tail DNA	1 µl
2 x Master Mix with enzymes and loading dye	12.5 µl
H ₂ O	9 µl
Total	25 µl

Competitor T's polymerase	
Primer set (10µM)	2.5 µl
Template: mouse tail DNA	1 µl
10 x Buffer	2.5 µl
MgCl (2.5 mM)	2 µl
dNTP (2.5 mM)	4 µl
PCR enzyme mixture	0.25 µl
H ₂ O	39.3 µl
Total	50 µl

PCR system:

T100 Thermal Cycler (BIO-RAD)

PCR cycling conditions:


Initial denaturation	95°C	3 min
Denaturation	95°C	15 sec
Annealing	55°C	15 sec
Extension	72°C	15 sec
Final Extension	72°C	5 min

35 cycles

Electrophoresis conditions:

Concentration of Agarose	1 %
Voltage	100 V
Time	30 min
Sample volume	10 µl

Results



The PCR specificity and easiness to use the kits were analyzed for FastGene® Optima^{HotStart ReadyMix} and a competitor. The results can be summarized as follows:

- FastGene® Optima^{HotStart ReadyMix} did not show non-specific amplicons.
- Competitor resulted in the presence of a large non-specific band.
- Amplification efficiency was equally high for both enzymes.
- PCR preparation was very easy with the ReadyMix.
- Electrophoresis was easily done since it already contained a loading dye.

FastGene® Optima for the Amplification of GC-Rich DNA

User's comment

"We saw less of non-specific amplification with GC-Rich fragments when using the FastGene® Optima. The sequencing results were at least comparable to those of competitor T's long amplicon reagent mixture."

The following data were provided by the courtesy of Ms. Ryoko Nakayama
Department of Pathology, Tsurumi University School of Dental Medicine, Japan.

Experimental conditions

The comparison of FastGene® Optima was done with a polymerase mixture designed for long amplifications from competitor T. Amplification of two GC-rich amplicons were done using the following reaction mixture:

PCR Reagents:

1. Competitor T enzyme Mix
2. FastGene® Optima (#LS28)

Template:

cDNA was created by reverse transcribing RNA isolated from a cell line of human origin

PCR Setup:

Competitor T's Polymerase mixture for long amplicons		
Primer set	1 µl	1 µl
Template: cDNA of a human cell line	2 µl	0.5 µl
2 x Buffer (GC-rich template)	12.5 µl	12.5 µl
dNTP	2 µl	2 µl
Enzyme mixture	0.2 µl	0.2 µl
H ₂ O	7.3 µl	8.8 µl
Total	25 µl	

FastGene® Optima		
Primer set	1 µl	1 µl
Template: cDNA of a human cell line	2 µl	0.5 µl
5 x Buffer (GC-rich template)	5 µl	5 µl
MgCl ₂	2.5 µl	2.5 µl
dNTP	1 µl	1 µl
FastGene® Optima	0.1 µl	0.1 µl
H ₂ O	13.4 µl	14.9 µl
Total	25 µl	

PCR product size and GC-content:

1. 1839 bp and 60.7%
2. 1260 bp and 64.3%


PCR cycling conditions:

Initial denaturation	94°C	5 min
Denaturation	94°C	30 sec
Annealing	58°C	30sec
Extension	72°C	2 min
Final Extension	72°C	7 min

35 cycles



Results



Competitor T's reagent and FastGene® Optima successfully amplified fragment (1) and (2). The specificity of both enzyme mixtures was high for both amplicons.

Amplified fragments were inserted into T-easy vector and transformed into E. coli. A selection of clones were send for sequencing to determine the number of mutations.

	Sequence	Clone No.	Determined sequence (bp)	Number of mutation	Number of mutation per Kb
1	Optima	18	1338	0	0.00
		20	422	1	2.37
		22	720	1	1.39
	Competitor T	6	1358	3	2.21
		8	608	4	6.58
2	Optima	1	1264	0	0.00
		1	1264	0	0.00
	Competitor T	11	381	1	2.62

Mutations were found in both enzyme mixtures. We analyzed however only a limited number of clones. Therefore, it is difficult to evaluate the accuracy. Nonetheless, FastGene® Optima had less mutations in when compared to competitor T.

Direct PCR from *Escherichia coli* colonies

User's comment

"The PCR products using the FastGene® Optima^{HotStart ReadyMix} with loading dye resulted in a clearer electrophoresis pattern when compared to competitor T's ready mix, since no smearing was observed. This smearing observed with the competitor's ready mix made the determination of the presence of the insert quite difficult. This is interesting, since the starting material (1 colony) was very similar. The starting amount was probably inadequate for the competitor's enzyme. For the FastGene® Optima^{HotStart ReadyMix} this was clearly not the case since one could clearly distinguish between colonies with insert (~800 bp) and without (~300 bp). Interestingly, good results were obtained using only 12.5 µl of reaction volume."

The following data were provided by an institutional customer in Japan.

Experimental conditions

PCR Reagents:

1. FastGene® Optima^{HotStart ReadyMix} with loading dye (#LS29)
2. Competitor T's enzyme ready mix with loading dye

Template:

24 transformed *Escherichia coli* colonies for each ready mix

PCR product size:

~800 bp (with insert) or
~300 bp (without insert)

PCR system:

Thermal Cycler Dice (Takara)

PCR Setup (for both polymerases):

FastGene® Optima ^{HotStart ReadyMix} with loading dye or Competitor T's ready mix with loading dye	
Primer set	0.5 µl (0.25µM)
Template: <i>E. coli</i> colony directly added to the reaction	1
2 x Master Mix with enzymes and loading dye	6.25 µl
H ₂ O	5.75 µl
Total	12,5 µl

PCR cycling conditions:

Initial denaturation	95°C	5 min
Denaturation	95°C	30 sec
Annealing	55°C	30 sec
Extension	72°C	2 min
Final Extension	72°C	5 min

} 35 cycles


Electrophoresis conditions:

Concentration of TAE- Agarose	1.2 %
Voltage	135 V
Time	30 min
Sample volume	5 µl




Results

FastGene® Optima^{HotStart ReadyMix} with loading dye



Competitor T's ready mix with loading dye



← 800 bp
← 300 bp

The electrophoresis pattern looked much better with FastGene® Optima^{HotStart ReadyMix}, when compared to the results obtained using the competitor's ready mix. Additionally, more primer dimers were observed when using the competitor's ready mix.

Amplification of complex DNA Templates

User's comment

„Almost no primer dimer were observed when using FastGene® Optima, especially when compared with competitor T's enzyme mixture. The amplification efficiency of the FastGene® Optima was high as well.“

The following data were provided by an institutional customer in Japan.

Experimental conditions

Comparison of the FastGene® Optima and a polymerase mixture of a competitor T was performed using the following conditions:

PCR Reagents:

1. Competitor T enzyme Mix
2. FastGene® Optima (#LS28)

Template:

cDNA was created by reverse transcribing RNA isolated from catshark liver tissue.

PCR system:

GeneAmp PCR System 9700 (ABI)

PCR Setup:

Competitor T's polymerase mixture for complex templates	
Primer set	1 µl (0.5µM)
Template: complex catshark's liver cDNA	0.5 µl (5ng)
10 x Buffer	5 µl
dNTP	4 µl
PCR enzyme mixture	0.2 µl
H ₂ O	39.3 µl
Total	50 µl

FastGene® Optima	
Primer set	1 µl (0.5µM)
Template: complex catshark's liver cDNA	0.5 µl (5ng)
5 x Buffer	10 µl
MgCl ₂	4 µl
dNTP	1.5 µl
FastGene® Optima	0.2 µl
H ₂ O	32.8 µl
Total	50 µl

PCR product size:

1030 bp

PCR cycling conditions:

Initial denaturation	94°C	2 min	} 35 cycles
Denaturation	94°C	1 min	
Annealing	59°C	45 sec	
Extension	72°C	1 min	
Final Extension	72°C	10 min	

Electrophoresis conditions:

Concentration of Agarose	1.2 %
Voltage	100 V
Time	35 min
Sample volume	20 µl

Results



Bulldog Bio

www.bulldog-bio.com
info@bulldog-bio.com

TEL 603-570-4248
FAX 603-766-0524

