

## BAMBANKER™ Serum-Free Cell Freezing Medium

For Research Use Only

Keep at 2~10°C

BAMBANKER™ is a serum-free medium for long-term freezing at -80° and preservation of valuable culture cells such as tumor cells and normal cells.

### Features

- Ready-to-use medium for preservation of cells,
- Usable without diluting,
- No program freezer is required
- Rapid and long-term freezing and preservation in a deep freezer (-80°)
- Serum-free

### Operating Procedure

- 1) Centrifuge the cells in logarithmic growth phase\*, and collect  $5 \times 10^5 \sim 1 \times 10^7$  cells.
- 2) Suspend the collected cells in 1 ml of this medium, and place the cells in cryo tubes for freezing and preservation. Then, freeze and preserve the cells at -80° without preliminary freezing. The cells frozen at -80° can be preserved subsequently in liquid nitrogen.
- 3) Thawing must be carried out quickly in a constant-temperature chamber or bath (37°C)
- 4) Bambanker contains DMSO, so if you like to remove this component please transfer 1 ml of your cell suspension to 10ml fresh culture media, mix well and remove the supernatant after centrifugation. Prepare and culture the cells as usual.

\*The cells must be frozen in logarithmic growth phase.

### Sterility test\*\*

- Endotoxin: Chromogenic substrate method
  - Mycoplasma: Fluorescent antibody method
  - Fungi and bacteria: As per Japanese Pharmacopoeia.
- (\*\*): Certificate of analysis will be issued upon request.)

### Applications (Cells preserved well by cryopreservation test):

P3U1 (mouse myeloma cell line), K562 (human leukemia cell line), human gastric epithelial cells, human  $\gamma\delta$ T cells.

Daudi (human B cell line), PC12 (rat-derived adrenal pheochromocytoma), human B cell line

OKT4 (mouse hybridoma), monkey B cell line, Activated lymphocyte derived from human peripheral blood, activated lymphocyte derived from mouse spleen

### Precautions

- ◆ Use only for research purposes, not for human use.
- ◆ Prior to using this medium, perform the confirmation test on the cells under study.
- ◆ The manufacturer is not responsible for any accident or damage caused by the use of this product.
- ◆ If you have any questions on the use of this medium, please contact the local distributor.

Catalog No.	Description	Package Size	Storage	Expiration period
BB01	<b>BAMBANKER™</b> Serum-Free Cell Freezing Medium, containing 10% DMSO	120mL	Keep at 2~10°C	2 years after manufacture
BB02		5 × 20 mL		
BB03		1 × 20 mL		

### Manufactured by:

**LYMPHOTEC Inc.** 5-26-9 Floral Building, Hakusan, Bunkyo-Ku, Tokyo 112-0001, Japan;

NIPPON GENETICS EUROPE GmbH; Binsfelder Strasse 77; 52351 Düren; HRB 4672

Bank data: UFJ Bank Limited: Identification Code: 30130700; Account number: 610762

Management board: Dr. Jürgen Lünzer; Tatsuji Hata, Kazuo Yamazaki

Value Added Tax Identification Number: DE 239977252