

결리분과  
리하리

Instruction manual

# MEGAquick-spin™ Total Fragment DNA Purification Kit

The Instruction Manual for fragment DNA Extraction from agarose gel, PCR product and enzymatic reaction using silica membrane.

REF	17286	Σ 50
REF	17287	Σ 200
REF	17288	Σ 200
RUO		
		
		Product info.
		



**INTRON**  
Biotechnology

[www.intronbio.com](http://www.intronbio.com)  
[www.intron-innoplex.com](http://www.intron-innoplex.com)

**INDEX**

Kit Information	Description	2
	Characteristics	2
	Kit Contents	3
	Storage	3
	Consideration Before Use	4
	Safety Information	4
	Additional Required Equipment	4
	Applications	5
	Quality Control	5
	Column Information	6
	Technical Assistance	6
Protocols	Protocol A (Gel DNA Extraction)	7
	Protocol B (PCR Purification, DNA Clean-up)	9
	Quick Guide	10
Additional information	Troubleshooting Guide	11
	Technical Advise	12
	Experimental Information	15
	Global Distributors	17

## DESCRIPTION

- The MEGAquick-spin™ Total Fragment DNA Purification Kit is designed to extract and purify DNA fragments of 60 bp ~ 20 kb from normal or low-melt agarose gels in either Tris acetate (TAE) or Tris borate (TBE), or to purify PCR products directly from a PCR amplification and DNA cleanup from other enzymatic reactions. Recovery is achieved up to 95%. PCR products are commonly purified to remove excess nucleotides and primers. This membrane-based system, which can bind up to 45 µg DNA, allows recovery of isolated DNA fragments or PCR products in as little as 20 minutes, depending on the number of samples processed and the protocol used. The purified DNA can be used for automated fluorescent DNA sequencing, cloning, labeling, restriction enzyme digestion or in vitro transcription/translation without further manipulation.
- The BNL Buffer are optimized for efficient recovery of DNA and removal of contaminants. As an added convenience from gel extraction procedures, the BNL Buffer contains a color indicator that allows easy monitoring of the solution pH for optimal DNA binding.
- The MEGAquick-spin™ Total Fragment DNA Purification Kit is ideal multi-functional (Gel extraction, PCR purification and DNA clean-up) product for isolation of fragment DNA.

## CHARACTERISTICS

- **Multifunction**  
Gel extraction, PCR purification and DNA clean-up
- **Improved recovery**  
Up to 95% recovery of ready-to-use DNA
- **Simple and easy process**  
Fast procedure, Cleanup of DNA up 60 bp ~ 20 Kb in three easy steps
- **Prevention of error**  
Easy determination of the optimal pH for DNA binding

**KIT CONTENTS**

Label	Description	Contents 50 Columns	Contents 200 Columns
<b>BNL Buffer<sup>1</sup></b>	Agarose gel lysis buffer	40 ml	140 ml
<b>Washing Buffer (concentrate)<sup>2</sup></b>	Washing buffer	10 ml	40 ml
<b>Elution Buffer</b>	Elution buffer	20 ml	20 ml
<b>MEGAquick-spin™ column (Blue column w/o Cap)</b>	Nucleic acid binding column	50 Col. (17286)	200 Col. (17287)
<b>MEGAquick-spin™ column (Clear Column with cap &amp; dark blue O-ring)</b>	Nucleic acid binding column		200 Col. (17288)
<b>Collection tube</b>	2 ml polypropylene tube	50 tubes	200 tubes



1 BNL Buffer contains chaotropic salts which are irritants. Take appropriate laboratory safety measures and wear gloves when handling.



2 Washing Buffer is supplied as concentrate. Add 40 ml (50 columns) or 160 ml (200 columns) per each bottles of ethanol (96~100%) according to the bottle label before use.

**STORAGE**

MEGAquick-spin™ Total Fragment DNA Purification Kit should be stored dry at room temperature (15–25°C). Under these conditions, MEGAquick-spin™ Total Fragment DNA Purification Kit can be stored for up to 24 months without showing any reduction in performance and quality. Check buffers for precipitate before use and redissolve at 37°C if necessary. The entire kit can be stored at 2–8°C, but in this case the buffers should be redissolved before use. Make sure that all buffers and spin columns are at room temperature when used. The term of validity is marked on the box.

### CONSIDERATION BEFORE USE

- A typical agarose gel slice is solubilized by adding 3 volumes of BNL Buffer to 1 volume of gel (e.g., 300 µl of BNL Buffer is added to 100 mg gel slice) and incubating at 55°C for 10 minutes. The high concentration of a chaotropic salt in BNL Buffer disrupts hydrogen bonding between sugars in the agarose polymer, allowing solubilization of the gel slice. In addition, the high salt concentration dissociates DNA binding proteins from the DNA fragments.
- If water is used for elution, make sure that its pH is between 7.0 and 8.5. Elution efficiency is dependent on pH and the maximum elution efficiency is achieved within this range. A pH <7.0 can decrease yield.

**Note :** Store DNA at -20°C when eluted with water, as DNA may degrade in the absence of a buffering agent.

### SAFETY INFORMATION

All chemicals should be considered as potentially hazardous. When working with chemicals, always wear a suitable lab coat and disposable glove. Some buffers contain the chaotropic salt which may be an irritant and carcinogen, so appropriate safety apparel such as gloves and eye protection should be worn. If a spill of the buffers occurs, clean with a suitable laboratory detergent and water.

If the liquid spill contains potentially infectious agents, clean the affected area first with laboratory detergent and water, then with a suitable laboratory disinfectant. Only persons trained in laboratory techniques and familiar with the principles of good laboratory practice should handle these products.



DO NOT add bleach or acidic solutions directly to the sample preparation waste.

### ADDITIONAL REQUIRED EQUIPMENT

- Agarose(iNtRON, Cta. No. 32032); scalpel
- Gel running buffer: TAE buffer or TBE buffer Electrophoresis Sterile
- Absolute ethanol
- Standard tabletop micro-centrifuge
- Micro-centrifuge tubes, sterile (1.5 ml)
- TE buffer (10 mM Tris-HCl, 0.1 mM EDTA; pH 8.0 - 8.5)

### APPLICATIONS

MEGAquick-spin™ Total Fragment DNA Purification Kit is designed for the efficient isolation of DNA fragments from TAE or TBE agarose gels or direct purification of PCR products. The purified DNA can be used for automated fluorescent DNA sequencing, cloning, restriction enzyme digestion and routinely performed DNA manipulation.

### QUALITY CONTROL

- In accordance with iNtRON's ISO-certified Total Quality Management System, each lot of MEGAquick-spin™ Total Fragment DNA Purification Kit is tested against predetermined specifications to ensure consistent product quality.
- The quality of the isolated fragment DNA was checked by agarose gel electrophoresis, and spectrophotometric determination.
- MEGAquick-spin™ column control  
The DNA binding capacity was tested by determining the recovery obtained with 20 µg of input fragment DNA. More than 70% recovery was obtained.
- Buffer control  
Conductivity and pH of buffers were tested and found to be within the pre-determined ranges

Buffer	Conductivity	pH
BNL	190 ~ 200 mS/cm	5.8 ~ 6.2
Washing B	14 ~ 16 mS/cm	7.4 ~ 7.8
Elution	450 ~ 700 µS/cm	8.0 ~ 8.5

### COLUMN INFORMATION

- The MEGAquick-spin™ Total Fragment DNA Purification Kit Spin Column

Column membrane <sup>1</sup>	Silica-based membrane
Spin Column <sup>1</sup>	Individually, in inserted in a 2.0 ml Collection Tube
Loading Volume	Maximum 800 µl
DNA Binding Capacity	Maximum 45 µg
Recovery	85 - 95% depending on the elution volume
Elution Volume	Generally, eluted with 30 – 200 µl of Elution Buffer



1. Do not store the Column packs under completely dried conditions. It may be affected to DNA binding capacity. The Spin Columns are stable for over 2 year under these conditions

### TECHNICAL ASSISTANCE

At iNtRON, we are proud of ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of iNtRON products. If you have any questions or experience any difficulties regarding the MEGAquick-spin™ Total Fragment DNA Purification Kit or iNtRON products in general, please do not hesitate to contact us.

**PROTOCOL A (Gel DNA Extraction)**

1. **Load and run the gel using an established protocol. DNA can be extracted from standard or low-melt agarose gels in TAE or TBE buffer.**

2. **After electrophoresis, cut out the interesting DNA fragment with a sharp scalpel or razor blade. Carefully take as much agarose gel as possible.**

**Note :** If sliced agarose gel put into BNL Buffer, the total volume may be increased. When highly concentrated BNL buffer is diluted, and it results low elution efficiency. Therefore, minimize the size of the gel slice by removing extra agarose.

**Note :** The gel slice may be stored at 4 °C or -20 °C for up to one week in a tightly closed tube under nuclease-free conditions before purification.

3. **Weigh the gel slice in a 1.5 ml tube. Add 3 volumes of BNL Buffer to 1 volume of gel (300 µl per 100 mg of agarose gel).**

**Note :** Add 300 µl of BNL Buffer to each 100 mg of gel. If more than 2% of agarose gel, add 6 volumes of BNL Buffer.

4. **Vortex the mixture and incubate at 55°C for 10 minutes or until the gel slice is completely dissolved. To help dissolving gel, vortex every 2 ~ 3 min during the incubation.**

**Note :** Vortex the tube every few minutes to increase the rate of agarose gel melting. Centrifuge the tube briefly at room temperature to ensure that the contents are at the bottom of the tube. Once the agarose gel is melted, the gel will not resolidify at room temperature.

**Note :** Completely solubilize agarose. For > 2% agarose gel, increase incubation time.



**Note :** If the color of the mixture is orange or violet, add 10 µl of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. The adsorption of DNA to the QIAquick membrane is efficient only at pH ≤7.5. BNL Buffer contains a pH indicator which is yellow at pH ≤7.5 and orange or violet at higher pH, allowing easy determination of the optimal pH for DNA binding.

5. **(Optional) For < 200 bp, add 1 gel volume of isopropanol to dissolved gel solution of the step 5 and mix well by pipetting several times. Do not centrifuge after mixing well.**



**Note :** For < 200 bp of DNA fragment, add 1 volume of isopropanol to 1 volume of gel, and mix well. If the agarose gel slice is 100 mg, add 100 µl of isopropanol.

When adding the isopropanol and mixing well by pipetting, small white pellet and clump should be formed. But never mind, and go to the following step. This step increases the yield of DNA fragment. For DNA fragment > 200 bp, adding isopropanol has no effect on yield.

6. **Place one MEGAquick-spin™ column in a Collection Tube for each dissolved gel mixture.**

- Transfer the dissolved gel mixture to the MEGAquick-spin™ column assembly.
- To bind DNA, apply the sample to the MEGAquick-spin™ column, and centrifuge for 1 min. Discard the flow-through after centrifuging and place the MEGAquick-spin™ column back in the same 2 ml collection tube.

**Note** : The maximum volume of the MEGAquick-spin™ column reservoir is 800 µl. For sample volumes of more than 800 µl, simply load and spin again.

- Add 700 µl of Washing Buffer to column and centrifuge at 13,000 rpm for 1 min. Discard the flow-through after centrifuging and place the MEGAquick-spin™ column back in the same 2 ml collection tube.

**Note** : If the DNA will be used for salt sensitive applications, such as blunt-end ligation and direct sequencing, repeat the step 4 using 500 µl of Washing Buffer.

- Centrifuge for 1 min at 13,000 rpm to dry the spin membrane.

**Note** : It is important to dry the spin membrane since residual ethanol may interfere with other reactions.

- Place the MEGAquick-spin™ column to a clean 1.5 ml microcentrifuge tube (not provided). Apply 30 ~ 100 µl of the Elution Buffer directly to the center of the column without touching the membrane with the pipette tip. Incubate at room temperature for 1 minute. Centrifuge for 1 minute at 13,000 rpm.

- Discard the MEGAquick-spin™ column and store the microcentrifuge tube containing the eluted DNA at -20°C.

**Note** : It is suggested to use at least 30 µl of the Elution Buffer to obtain best result.

**Note** : Ensure that the elution buffer is dispensed directly onto the MEGAquick-spin™ membrane for complete elution of bound DNA. The average eluate volume is 48 µl from 50 µl elution buffer volume, and 28 µl from 30 µl Elution Buffer.

**Note** : Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5. When using water, make sure that the pH value is within this range, and store DNA at -20°C as DNA may degrade in the absence of a buffering agent. The purified DNA can also be eluted in TE buffer (10 mM Tris·Cl, 1 mM EDTA, pH 8.0), but the EDTA may inhibit subsequent enzymatic reactions.

**PROTOCOL B (PCR Purification, DNA Clean-up)**

1. **Amplify target sample using standard amplification conditions. (or Prepare the DNA mixture enzymatically reacted for clean-up)**
2. **Add an 5 volume of BNL Buffer to the PCR reaction product, and mix well by vortexing. If the PCR product is 20  $\mu$ l, add 100  $\mu$ l of BNL buffer to the PCR tube directly.**

**Note :** Centrifuge the tube briefly at room temperature to ensure that the contents are at the bottom of the tube.



**Note :** If the color of the mixture is orange or violet, add 10  $\mu$ l of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. The adsorption of DNA to the MEGAquick-spin™ membrane is efficient only at pH  $\leq$ 7.5. BNL Buffer contains a pH indicator which is yellow at pH  $\leq$ 7.5 and orange or violet at higher pH, allowing easy determination of the optimal pH for DNA binding.

3. **(Optional) For < 200 bp, Add 1.5 volume of isopropanol to the sample and mix well by pipetting several times. Do not centrifuge after mixing well.**



**Note :** For < 200 bp, Add 1.5 volume of isopropanol, and mix well. If the PCR product is 20  $\mu$ l, add 100  $\mu$ l of BNL Buffer and 150  $\mu$ l of isopropanol. This step increases the yield of DNA fragment.

4. **Place one MEGAquick-spin™ column in a Collection Tube for each DNA gel mixture.**
5. **Transfer the DNA mixture to the MEGAquick-spin™ column assembly.**
6. **To bind DNA, apply the sample to the MEGAquick-spin™ column, and centrifuge for 1 min. Discard the flow-through after centrifuging and place the MEGAquick-spin™ column back in the same 2 ml collection tube.**  
**Note :** The maximum volume of the MEGAquick-spin™ column reservoir is 800  $\mu$ l. For sample volumes of more than 800  $\mu$ l, simply load and spin again.
7. **Add 700  $\mu$ l of Washing Buffer to column and centrifuge at 13,000 rpm for 1 min. Discard the flow-through after centrifuging and place the MEGAquick-spin™ column back in the same 2 ml collection tube.**  
**Note :** If the DNA will be used for salt sensitive applications, such as blunt-end ligation and direct sequencing, repeat the step 4 using 500  $\mu$ l of Washing buffer.
8. **Centrifuge for 1 min at 13,000 rpm to dry the spin membrane.**

**Note :** It is important to dry the spin membrane since residual ethanol may interfere with other reactions.


- Place the MEGAquick-spin™ column to a clean 1.5 ml microcentrifuge tube (not provided). Apply 30 ~ 100 µl of the Elution Buffer directly to the center of the column without touching the membrane with the pipette tip. Incubate at room temperature for 1 minute. Centrifuge for 1 minute at 13,000 rpm.
- Discard the MEGAquick-spin™ column and store the microcentrifuge tube containing the eluted DNA at -20°C.

**Note :** It is suggested to use at least 30 µl of the Elution Buffer to obtain best result.

**Note :** Ensure that the elution buffer is dispensed directly onto the MEGAquick-spin™ membrane for complete elution of bound DNA. The average eluate volume is 48 µl from 50 µl elution buffer volume, and 28 µl from 30 µl elution buffer.

**Note :** Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5. When using water, make sure that the pH value is within this range, and store DNA at -20°C as DNA may degrade in the absence of a buffering agent. The purified DNA can also be eluted in TE buffer (10 mM Tris-Cl, 1 mM EDTA, pH 8.0), but the EDTA may inhibit subsequent enzymatic reactions.

### Quick Guide



## TROUBLESHOOTING GUIDE

Problem	Possible Cause	Recommendation
Low or no yield	Washing buffer did not contain ethanol	<ul style="list-style-type: none"> <li>Ethanol must be added to Washing buffer before use.</li> </ul>
	Inappropriate elution buffer	<ul style="list-style-type: none"> <li>DNA will only be eluted in low salt buffer or water.</li> </ul>
	Incorrect volume of BNL buffer	<ul style="list-style-type: none"> <li>Verify that a correct volume of BNL buffer was added to the gel slice.</li> </ul>
	Gel slice incompletely solubilized	<ul style="list-style-type: none"> <li>After addition of BNL buffer to the slice, mix by vortexing the tube every 2 minutes during the 55 °C incubation.</li> </ul>
DNA does not perform well, e.g., in enzyme reaction, ligation, sequencing reactions	Cloudy and gelatinous appearance of sample mixture after addition of isopropanol	<ul style="list-style-type: none"> <li>This may be due to salt contamination, and will be disappeared by mixing the sample. Alternatively, the gel slice may not be completely solubilized. The concentration of gel may be above 2%. In this case, apply the 6 volume of BNL buffer to gel slice, and melt the gel completely.</li> </ul>
	Salt concentration in eluate too high	<ul style="list-style-type: none"> <li>Modify the washing step by incubating the column for 5min at RT after adding 700 µl of Washing Buffer and the centrifuge.</li> </ul>
	Eluate contaminated with agarose	<ul style="list-style-type: none"> <li>If The gel slice is incompletely solubilized or overweighed, increased the incubation time.</li> </ul>
	Eluate contains residual ethanol	<ul style="list-style-type: none"> <li>Ensure that the wash flow-through is drained from the collection tube and that the column is then centrifuged at 13,000 rpm for 1min.</li> </ul>
BNL Buffer become violet color.	Eluate contain primer-dimers	<ul style="list-style-type: none"> <li>Primer dimers formed are longer than 50bp, and are not completely removed.</li> <li>After binding step, wash the column with 750 µl of a 35% guanidine hydrochloride aqueous solution. Follow with the washing, and elution step as in the protocols.</li> </ul>
	pH of electrophoresis buffer too high	<ul style="list-style-type: none"> <li>The electrophoresis buffer has been repeatedly used or incorrectly prepared, resulting in a sample pH that exceeds the buffering capacity of BNL Buffer and leads to inefficient DNA binding. Add 0.1 volume of 3M sodium acetate, pH 5.0, to the sample and mix.</li> </ul>

## TECHNICAL ADVICE

### ◆ Principle of gel analysis

Gels allow separation and identification of nucleic acids based on charge migration. Migration of nucleic acid molecules in an electric field is determined by size and conformation, allowing nucleic acids of different sizes to be separated. However, the relationship between the fragment size and rate of migration is non-linear, since larger fragments have greater frictional drag and are less efficient at migrating through the polymer.

Agarose gel analysis is the most commonly used method for analyzing DNA fragments between 0.1 and 25 Kb. Other specialized analytical gel methods exist for analyzing extremely large or small DNA molecules. Detailed information on all types of analytical gels can be found in current molecular biology manuals.

### ◆ Agarose concentration

The concentration of agarose used for the gel depends primarily on the size of the DNA fragments to be analyzed. Low agarose concentrations are used to separate large DNA fragments, while high agarose concentrations allow resolution of small DNA fragments

#### Concentration of agarose used for separating DNA of different sizes

Agarose conc. (% w/v)	DNA fragment range (Kb)	Agarose conc. (% w/v)	DNA fragment range (Kb)
0.3	5 – 60	0.5	1 – 30
0.7	0.8 – 12	1.0	0.5 – 10
1.2	0.4 – 7	1.5	0.2 – 3
2.0	0.05 – 2		

### ◆ Electrophoresis buffers


The most commonly used buffers for agarose gel electrophoresis are TBE (Tris-borate–EDTA) and TAE (Tris-acetate–EDTA). Although more frequently used, TAE has a lower buffering capacity than TBE and is more easily exhausted during extended electrophoresis. TBE gives better resolution and sharper bands, and is particularly recommended for analyzing fragments <1 Kb.

The drawback of TBE is that the borate ions in the buffer form complexes with the cis-diol groups of sugar monomers and polymers, making it difficult to extract DNA fragments from TBE gels using traditional methods.

### ◆ Effect of pH on DNA Binding

The MEGAquick-spin™ column is uniquely adapted to purify DNA from both aqueous solutions and agarose gels, and up to 45 µg DNA can bind to each MEGAquick-spin™ column. The binding buffers in MEGAquick-spin™ Total Fragment DNA Purification Kit provide the correct salt concentration and pH for adsorption of DNA to the MEGAquick-spin™ column. The adsorption of nucleic acids to silica surfaces occurs only in the presence of a high concentration of chaotropic salts, which modify the structure of water. Adsorption of DNA to silica also depends on pH. Adsorption is typically 95% if the pH is  $\leq 7.5$ , and is reduced drastically at higher pH (Figure 1). If the loading mixture pH is  $>7.5$ , the optimal pH for DNA binding can be obtained by adding a small volume of 3 M sodium acetate, pH 5.0.

BNL Buffer contains an integrated pH indicator allowing easy determination of the optimal pH for DNA binding. DNA adsorption requires a pH 7.5, and the pH indicator in the buffers will appear yellow in this range. If the pH is  $>7.5$ , which can occur if during agarose gel electrophoresis, the electrophoresis buffer had been used repeatedly or incorrectly prepared, or if the buffer used in an enzymatic reaction is strongly basic and has a high buffering capacity, the binding mixture turns orange or violet (Fig. 1). This means that the pH of the sample exceeds the buffering capacity of BNL Buffer and DNA adsorption will be inefficient. In these cases, the pH of the binding mixture can easily be corrected by addition of a small volume of 3 M sodium acetate, pH 5.0, before proceeding with the protocol.



**Fig. 1. pH dependence of DNA adsorption MEGAquick-spin™ column.**

1 µg of a 2.9 Kb DNA fragment was adsorbed at different pHs and eluted with Elution Buffer. The graph shows the percentage of DNA recovery, reflecting the relative adsorption efficiency, versus pH of adsorption. The pH indicator dye in MEGAquick-spin™ Total Fragment DNA Purification Kit identifies optimal pH for DNA binding.

### ◆ Low Salt Elution

Elution efficiency is strongly dependent on the salt concentration and pH of the Elution Buffer. Contrary to adsorption, elution is the most efficient under basic conditions and low salt concentrations. DNA is eluted with 50 or 30  $\mu\text{l}$  of the provided Elution Buffer, or water. The maximum elution efficiency is achieved between pH 7.0 and 8.5. When using water to elute, make sure that the pH is within this range. In addition, DNA must be stored at  $-20^{\circ}\text{C}$  when eluted with water since DNA may degrade in the absence of a buffering agent. Elution with TE buffer (10 mM Tris-Cl, 1 mM EDTA, pH 8.0) is possible, but not recommended because EDTA may inhibit subsequent enzymatic reactions.


### ◆ DNA yield and concentration

DNA yield depends on the following three factors: the volume of elution buffer, how the buffer is applied to the column, and the incubation time of the buffer on the column. 100–200  $\mu\text{l}$  of Elution Buffer completely covers the MEGAquick-spin™ column, ensuring maximum yield, even when not applied directly to the center of the membrane. Elution with  $\leq 50$   $\mu\text{l}$  requires the buffer to be added directly to the center of the membrane, and if elution is done with the minimum recommended volume of 30  $\mu\text{l}$ , an additional 1 minute incubation is required for optimal yield. DNA will be up to 1.7 times more concentrated if the MEGAquick-spin™ column is incubated for 1 minute with 30  $\mu\text{l}$  of Elution Buffer, than if it is eluted in 50  $\mu\text{l}$  without incubation

## EXPERIMENTAL INFORMATION

### ◆ Yields of various sizes of Fragment DNA

MEGAquick-spin™ Total Fragment DNA Purification Kit shows improved DNA recovery from short size of DNA to long length fragment (20 Kb).




**Fig. 2. Yield of fragment DNA**

- Fragment DNA Size : 87 bp ~ 20 Kb,

Lane M1 & M2, DNA marker; lane 1, Before purification; lane 2, after purification with older Product; lane 3, after purification with MEGAquick-spin™ Total Fragment DNA Purification Kit.


### ◆ Comparative Test of DNA Recovery

MEGAquick-spin™ Total Fragment DNA Purification Kit shows the performance of competitive advantage of DNA extraction from all fragment size and amount.



**Fig. 3. Comparative performance test of fragment DNA purification (Gel Ext.)**


- Fragment DNA Size : 108 bp ~ 20 Kb,



**Fig. 4. Comparative performance test of DNA purification (PCR purification)**

- Fragment DNA Size : 1kb fragment (linearized PCR fragment)
- Start DNA amount : 100 ng ~ 2000 ng
- Average recovery : MEGAquick-spin™ Total - 85 ~ 97%, Supplier A - 75 ~ 85%, Supplier B - 80 ~ 95%

#### ◆ Suitable for Down-stream Operations



**Fig. 5. Reliable Long Read Lengths in Sequencing**

High quality sequencing data of 4.5 Kb length DNA fragment purified with iNtRON's MEGAquick-spin™ Total Fragment DNA Purification Kit .

## ◆ Molecular Reagent

### Australia

Scientifix Pty Ltd.  
Phone : +61 3 85405900  
Fax : +61 3 9548 7177  
URL : <http://www.scientifix.com.au>

### Belgium

European Biotech Network  
Phone : +32 4 3884398  
Fax : +32 4 3884398  
URL : <http://www.euro-bio-net.com>

### Canada

FroggaBio  
Phone : +1 416 736 8325  
Fax : +1 416 736 3399  
URL : <http://www.froggabio.com>

### China

Chinagen Inc.  
Phone : +86 (0)755 26014525  
Fax : +86 (0)755 26014527  
URL : <http://www.chinagen.com.cn>

### China - Hong Kong

Tech Dragon Limited  
Phone : +852 2646 5368  
Fax : +852 2646 5037  
URL :  
<http://www.techdragon.com.hk>

### Egypt

Biovision Egypt Co.  
Phone : +20 119007908  
Fax : +20 223204509  
Email : [biovision.egypt@gmail.com](mailto:biovision.egypt@gmail.com)

### France

EUROMEDEX  
Phone : +33 3 88 18 07 22  
Fax : +33 3 88 18 07 25  
URL : <http://www.euromedex.com>

### Germany

HISS Diagnostics GmbH.  
Phone : +49 761 389 490

Fax : +49 761 202 0066  
URL : <http://www.hiss-dx.de>

### Hungary

Bio-Kasztel Kft.  
Phone : +36 1 381 0694  
Fax : +36 1 381 0695  
URL : <http://www.kasztel.com>

### India

Biogene  
Phone: +91 11  
42581008/25920048  
fax: +91 11 42581260  
URL : <http://www.biogene-india.com>

### Indonesia

CV.Kristalindo Biolab  
Phone : +62 31 5998626  
Fax : +62 31 5998627  
Email : [kutama@indo.net.id](mailto:kutama@indo.net.id)

### Iran

NANOMEHR CO.  
Phone : +98 21 4432 3682  
Fax : +98 21 4432 3684  
URL : <http://www.nanomehr.ir>

### Israel

Talron Biotech Ltd.  
Phone : +972 8 9472563  
Fax : +972 8 9471156  
URL : <http://www.talron.co.il>

### Italy

Li StarFISH S.r.l  
Phone : +39 02 92150794  
Fax : +39 02 92157285  
URL : <http://www.listarfish.it>

### Japan

Cosmo Bio Co.,LTD.  
Phone : +81 3 5632 9617  
Fax : +81 3 5632 9618  
URL :  
<http://www.cosmobio.co.jp>

### Netherlands

Goffin Molecular Technologies B.V.  
Phone : +31 76 508 6000  
Fax : +31 76 508 6086  
URL : <http://www.goffinmeyvis.com>

### New Zealand

Ngaiio Diagnostics Ltd  
Phone : +64 3 548 4727  
Fax : +64 3 548 4729  
URL : <http://www.ngaiio.co.nz>

### Spain

LABOTAQ, S.C  
Phone : +34 954 31 7216  
Fax : +34 954 31 7360  
URL : <http://www.labotaq.com>

### Taiwan

Asian Life Science Co. Ltd.  
Phone : +886 2 2998 6239  
Fax : +886 2 8992 0985  
URL : <http://www.asiansci.com.tw>

### Taiwan

Hong-jing Co., Ltd.  
Phone : +886 2 3233 8585  
Fax : +886 2 3233 8686  
URL :  
<http://www.hongjing.com.tw>

### Thailand

Pacific Science Co. Ltd.  
Phone : +66 2 433 0068  
Fax : +66 2 434 2609  
URL :  
<http://www.Pacificscience.co.th>

### Turkey

BIOCEM Ltd. Co.  
Phone : +90 212 534 0103  
Fax : +90 212 631 2061  
URL : <http://www.biocem.com.tr>



### ◆ Molecular Reagent

#### United Kingdom

CHEMBIO LTD.  
Phone : +44 208 123 3116  
Fax : +44 800 007 3116  
URL : <http://www.chembio.co.uk>

#### U.S.A.

Boca Scientific  
Phone : +1 561 995 5017  
Fax : +1 561 995 5018  
URL : <http://www.bocascientific.com>

#### Vietnam

VIETLAB Co., Ltd  
Phone : +844 37821739  
Fax : +844 37821738  
Email : [info@vietlab.vn](mailto:info@vietlab.vn)

### ◆ Molecular Diagnosis

#### Iran

Sina Bio Medical Chemistry Co.  
Phone : +98 21 2244 2488  
Fax : +98 21 2244 0888  
URL : <http://www.sinabiomedical.com>

#### Kazakhstan

BioHim Pribor  
Phone : +7 727 278 23 16  
Fax : +7 727 269 2791  
Email : [biohimpribor@mail.ru](mailto:biohimpribor@mail.ru)

#### Spain

EUROVET VETERINARIA S.L.  
Phone : +34 91 8841374  
Fax : +34 918875465  
URL : <http://www.euroveterinaria.com>

#### Iran

Sina Bio Medical Chemistry Co.  
Phone : +98 21 2244 2488  
Fax : +98 21 2244 0888  
URL : <http://www.sinabiomedical.com>

#### Kazakhstan

BioHim Pribor  
Phone : +7 727 278 23 16  
Fax : +7 727 269 2791  
Email : [biohimpribor@mail.ru](mailto:biohimpribor@mail.ru)

#### Spain

EUROVET VETERINARIA S.L.  
Phone : +34 91 8841374  
Fax : +34 918875465  
URL : <http://www.euroveterinaria.com>

### ◆ Molecular Reagent / Molecular diagnosis

#### Austria

Anopoli Biomedical Systems  
Phone : +43 2773 42564  
Fax : +43 2773 44393  
URL : <http://www.anopoli.com>

#### Jordan / Iraq

Genetics Company  
Phone : +962 6 5536402  
Fax : +962 6 5536398  
URL : <http://www.genetics-jo.com>

#### Malaysia

NHK BIOSCIENCE SOLUTIONS  
SDN  
Phone : +60 3 7987 8218  
Fax : +60 3 7987 8213  
URL : <http://www.nhkbioscience.com>

#### Mongolia

SX Biotech Co., Ltd.  
Phone : +976 5006 0677  
Fax : +976 7011 1767  
Email : [zanaa@sxbiotech.mn](mailto:zanaa@sxbiotech.mn)

#### Pakistan

HR BIO SCIENCES  
Phone : +92 42 37247650  
Fax : +92 42 37247650  
Email : [hribiosciences@yahoo.com](mailto:hribiosciences@yahoo.com)

#### Philippines

Hebborn Analytics INC.  
Phone : +632 461 7173  
Fax : +632 418 5877  
Email : [hebborn@pltdsl.com](mailto:hebborn@pltdsl.com)

#### Romania

S.C. Bio Zyme S.R.L.  
Phone : +40 264 52 32 81  
Fax : +40 264 52 32 81  
URL : <http://www.biozyme.ro>

#### Switzerland

LucernaChem AG  
Phone : +41 (0)41 420 9636  
Fax : +41 (0)41 420 9656  
URL : <http://www.lucerna-chem.ch>

#### Tunisia

RIBO Pharmaceutique &  
Diagnostique  
Phone : +216 71981095  
Fax : +216 71981473  
Email : [nbopharma@topnet.tn](mailto:nbopharma@topnet.tn)

#### U.S.A.

Bulldog Bio Inc.  
Phone : +1 603 570 4248  
Fax : +1 603 766 0524  
URL : <http://www.bulldog-bio.com>

## Customer & Technical Service

Do not hesitate to ask us any question

### Contact to us



shop.intronbio.com

Tel : +82-505-550-5600

Fax : +82-505-550-5660

Mail : intronbio@intronbio.com

### Near your partner

You can find your partners, iNTRON Distributor in Page.

### Multifunction

Gel extraction, PCR purification and DNA clean-up

### Improved recovery

Up to 95% recovery of ready-to-use DNA

### Fast procedure

Cleanup of DNA up 60 bp ~ 20 Kb in three easy steps

### Prevention of error

Easy determination of the optimal pH for DNA binding



Rm 701~ 704, Jung-Ang Induspia B/D  
Sangdaewon-dong, Joongwon-gu,  
Sungnam-si, Kyeonggi-do, Korea  
iNTRON Biotechnology, Inc.

Trademarks: iNTRON, DNA-spin™, DNA-midi™, DNA-maxi™, PCRquick-spin™, MEGA-spin™, MEGAquick-spin™, MEGA-bead™, PROBER™, G-DEX™, G-spin™, Viral Gene-spin™, easy-spin™, RNA-spin™, easy-BLUE™, easy-RED™, WEST-one™, WEST-ZOL™, PRO-PREP™, SMART™, PRO-MEASURE™, Genelator™, F-Detector™, Broad-Way™, PRO-STAIN™, pLUG, Maxime™, i-Taq™, i-StarTaq™, i-MAX™, i-StarMAX™, RedSafe™, Muta-Direct™, e-Mycos™, M-Solution™, CENDORI™, VeTeK™, iNNOPLEX™, GxN™, teleFAXgene™, CLP™ and IQeasy™ is a trademark of iNTRON Biotechnology, Inc.

iNTRON kits are intended for research use only. Prior to using them for other purposes, the user must validate the system in compliance with the applicable law, directives, and regulations.

The PCR process is covered by patents issued and applicable in certain countries. iNTRON Biotechnology, Inc. does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.

© 2011 iNTRON, all rights reserved.