



Application

Comparative evaluation of FastGene® Scriptase II and other company's kit by RT-PCR

Category

Reverse transcriptase

Product

FastGene® Scriptase II (LS53)

Manufacturer

NIPPON Genetics EUROPE

The following data has been published due to the kindness of Mr. Tamura, a member of the Japanese public research institute.

Overview



FastGene® Scriptase II (LS53)

In the order to evaluate the performance of FastGene® Scriptase II, reverse transcription was performed on the RNA extracted from HeLa cells with this product and the other company's reverse transcriptase. The obtained cDNA was amplified by PCR for the target gene and confirmed by electrophoresis. The performance of FastGene® Scriptase II was evaluated by comparing the band signals.

Experiment

• RNA extraction

Sample: HeLa cells were cultured on 35 mm dish and used at 80-90% confluent time point.

RNA extraction: Ribozol (MS Technosystems) (N580-100ML)

Absorbance measurement: U-3900 (HITACHI) (U-3900)

Result: A260/280 = 1.72, yield = 41.5 µg

• Reverse transcriptase reaction

Total RNA was diluted to 1, 10, 100 ng/µL. 1 µL of each dilution was added as Template RNA. ※ 0 ng is 1 µL H₂O

• Primer

Gene specific primer (GAPDH) : CTCTTCCTCTGTGCTCTTGC

FastGene® Scriptase II

Template RNA (0, 1, 10, 100 ng/µL)	: 1 µL
Gene specific primer (GAPDH)	: 2 µL (2 pmol)
dNTP Mixture (2 mM each)	: 2 µL
Sterile water (RNase free)	: 7.5 µL
total	: 12.5 µL
↓	
65 °C	5 min
↓	
4 °C	5 min
↓	
add	
5X FastGene® Scriptase II buffer	: 4 µL
0.1 M DTT (included in kit)	: 2 µL
FastGene® Scriptase II	: 1 µL
↓	
42 °C	50 min
↓	
70 °C	15 min
↓	
keep at 22 °C	

Reference: Template total RNA recommended amount: 1 ng-5 µg

company T, RT reaction kit

Template RNA (0, 1, 10, 100 ng/µL)	: 1 µL
Gene specific primer (GAPDH)	: 2 µL (2 pmol)
2 mM dNTP	: 5 µL
RNase free dH ₂ O	: 2 µL
total	: 10 µL
↓	
65 °C	5 min
↓	
4 °C	5 min
↓	
add	
5X buffer	: 4 µL
Reverse transcriptase	: 1 µL
RNase free dH ₂ O	: 5 µL
↓	
42 °C	50 min
↓	
70 °C	15 min
↓	
keep at 22 °C	

Reference: Template total RNA recommended amount: 5 µg or less



PCR conditions

● Reaction composition

10 μM Forward Primer	: 2.5 μL
10 μM Reverse Primer	: 2.5 μL
2 mM dNTPs (TOYOBO)	: 5 μL
5X Phusion HF Buffer	: 10 μL
cDNA	: 2 μL
sterile water	: 27.5 μL
Phusion DNA Polymerase (NEB)	: 0.5 μL
total	: 50 μL

● Program

95 °C	2 min	
↓		
95 °C	10 sec	} 30 cycles
65 °C	30 sec	
72 °C	30 sec	
↓		
72 °C	3 min	
↓		
hold	22 °C	

● Primers for-GAPDH

Forward Primer : CCACAGTCCATGCCATCAC
 Reverse Primer : CCATGAGGTCCACCACCC
 Amplification product size: 500bp

Electrophoresis condition

Agarose gel: 1% agarose gel


Sample: PCR reaction solution (15 μL)

Running buffer: 1xTAE

Migration conditions: 100 V, 30 min

Marker: 7 μL (Gene Ladder Wide 1) (Nippon Gene) (313-06961)

Result



We were able to detect PCR bands by using Scriptase II and a small amount of total RNA.



Customer's comment

In order to analyze the mRNA expression level of cultured cells, I was searching for a reverse transcriptase that is convenient and can analyze a large amount of samples.

I am satisfied with reverse transcription with a small amount of total RNA.