



TaqDog Premium DNA Polymerase

Products

CAT. #	Description	Volume
TQDG010	TaqDog Premium DNA Polymerase, with 10x Reaction Buffer (no dNTP), 1000U	1,000 Units (200 µL)
TQDG050	TaqDog Premium DNA Polymerase, with 10x Reaction Buffer (no dNTP), 5,000U	5,000 Units (1.0 mL)
TQDG100	TaqDog Premium DNA Polymerase, with 10x Reaction Buffer (no dNTP), 10,000U	10,000 Units (2 x 1.0 mL)
TQDG500	TaqDog Premium DNA Polymerase, with 10x Reaction Buffer (no dNTP), 50,000U	50,000 Units (10 x 1.0 mL)

Product Information

TaqDog Premium DNA Polymerase generates a consistent PCR product across a range of DNA templates, including mouse or human, plasmid and bisulfite converted DNA. The Taq polymerase leaves an adenine overhang at the five-prime end making the PCR product ideal for TA cloning.

Extension: 3,000 base-pair extension has been confirmed.

Speed: Extension speeds up to 1 Kbp / minute.

Product Components

	50,000 U	10,000 U	5000 U	1,000 U
Taq Polymerase	10 x 1 ml	2 x 1 ml	1 x 1 ml	1 x 200 µl
10X PCR Buffer	5 x 50 ml	1 x 50 ml	1 x 50 ml	4 x 10 ml
MGSO ₄	2 x 50 ml	1 x 50 ml	1 x 10 ml	2 x 1.5 ml
GC - Extender	3 x 50 ml	3 x 10 ml	2 x 10 ml	2 x 1.5 ml

Taq Storage Buffer: (50 mM Tris-HCl (pH 8.0), 100 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 1 % Triton X-100, 50 % glycerol)

10X PCR Buffer: 100 mM Tris-HCl (pH 8.8), 100 mM KCl, 100 mM (NH₄)₂SO₄, 1% Triton X-100

All components are to be stored at -20 °C



Protocol

1. Add following components to a PCR tube on ice:

Note: DNA amount, primer amount and GC Extender use/volume should be optimized.

Component	50 µl reaction
Template DNA (up to 500 ng)	up to 10 µl
10 µM Forward Primer	0.5 µl
10 µM Reverse Primer	0.5 µl
10X PCR Buffer	5 µl
50 mM MgSO ₄	1.5 µl
10 mM dNTP mix	1 µl
Taq DNA Polymerase	0.2 µl
GC-Extender (Optional)*	5 µl – 15 µl
Nuclease-free H ₂ O	up to 50 µl

*GC-Extender is intended for difficult to amplify templates with GC > 60 %, or for templates longer than 2.5 Kbp. Volume to use (between 5 µl – 15 µl) should be optimized for each template and primer combination.

2. Mix contents (either through pipetting multiple times or briefly vortexing).

3. Centrifuge briefly and proceed to PCR amplification steps.

4. PCR amplification:

Note: Number of cycles, Annealing temperature and extension should be optimized for each template and primer set.

Stage	Step	Temperature (°C)	Time	Cycles
1	Denature	94	10 sec - 30 sec	20-30
	Anneal	55 - 60	10 sec - 30 sec	
	Extend	72	30 sec	
2	Final extension step	72	7 min	1
3	Holding step	4	∞	1

5. Product is ready for downstream analysis. For long-term storage store at -20C.

END OF PROTOCOL