

# C-Chip

The World's First True  
Disposable Hemocytometer



Precise volume control

Easy to use



# C-Chip is



a disposable plastic hemocytometer, which is a significant improvement from the conventional glass hemocytometer.

- Guarantees repeatability and reliability through complete quality control
- Two exact 100  $\mu\text{m}$  height chambers do not need coverslips
- Eliminates the need of tedious washing
- Reduces exposure to infectious samples and contamination
- Sturdy and strong as it is made from very right durable plastic


## Applications

**Blood analysis (Hematology)** | Blood cell counting

**Cell Culture** | Cell concentration measurement / Cell viability test

**Microbiology** | Bacteria & Fungal spore counting

**IVF, IUI** | Sperm counting



No.	Grid	Depth
DHC-N01	Neubauer Improved	100 $\mu\text{m}$

## Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm

Chamber depth : 100, 200, 10  $\mu\text{m}$

Chamber volume : 10  $\mu\text{L}$

