


BioReady™ rTaq DNA Polymerase

PCR Enzyme. **3¢/unit** What More Could You Ask For?

Amplify Your Results

Cost effective – BioReady Taq is a great choice for your everyday PCR needs. This recombinant DNA polymerase is optimized for improved efficiency and can amplify between 50 bp and 5,000 bp. When used with any of our thermal cyclers, protocols can be completed in less than 1 hour. BioReady Taq is manufactured under stringent GMP control, and though we don't recommend it, this ultra-pure enzyme extract is still fully active after days at 37°C.



As low as
3-cents/unit!

SPECIFICATIONS

Species of Origin	Thermus auaticus
Enzyme Type	Thermostable DNA polymerase
Form	Recombinant - Engineered
5' – 3' Exonuclease Activity	Yes
"A" Overhang	Yes
3' – 5' Exonuclease Activity	No
Reaction Buffer	10x, 15mM MgCl ₂
Storage	1 year at -20°C
Concentration	5 units/μl
DNase Free	Yes
RNase Free	Yes
Protease Free	Yes

ITEM #	DESCRIPTION	PRICE
BSA12L050	BioReady rTaq, 5,000 units, 3-cents/unit	\$150
BSA12L020	BioReady rTaq, 2,000 units, 4-cents/unit	\$80
BSA12L010	BioReady rTaq, 1,000 units, 5-cents/unit	\$50