


BioReady™ rTaq DNA Polymerase

PCR Enzyme. **3¢/unit** What More Could You Ask For?

Amplify Your Results

Cost effective – BioReady Taq is a great choice for your everyday PCR needs. This recombinant DNA polymerase is optimized for improved efficiency and can amplify between 50 bp and 5,000 bp. When used with any of our thermal cyclers, protocols can be completed in less than 1 hour. BioReady Taq is manufactured under stringent GMP control, and though we don't recommend it, this ultra-pure enzyme extract is still fully active after days at 37°C.



As low as
3-cents/unit!

SPECIFICATIONS

| | |
|------------------------------|-----------------------------|
| Species of Origin | Thermus auaticus |
| Enzyme Type | Thermostable DNA polymerase |
| Form | Recombinant - Engineered |
| 5' – 3' Exonuclease Activity | Yes |
| "A" Overhang | Yes |
| 3' – 5' Exonuclease Activity | No |
| Reaction Buffer | 10x, 15mM MgCl ₂ |
| Storage | 1 year at -20°C |
| Concentration | 5 units/μl |
| DNase Free | Yes |
| RNase Free | Yes |
| Protease Free | Yes |

| ITEM # | DESCRIPTION | PRICE |
|-----------|--|-------|
| BSA12L050 | BioReady rTaq, 5,000 units, 3-cents/unit | \$150 |
| BSA12L020 | BioReady rTaq, 2,000 units, 4-cents/unit | \$80 |
| BSA12L010 | BioReady rTaq, 1,000 units, 5-cents/unit | \$50 |