


Cell Fusion Electrodes

Electrodes for Hybridoma Production, Mammalian Cell Fusion, etc.


MS Stand Model Chamber Platinum Electrodes & Container



[CUY497P2](#) [CUY497P3](#) [CUY497P4](#) [CUY499L](#)

- [CUY497P2](#) MS Stand Model Chamber Platinum Electrode
2mm gap L80xW2xH5mm 0.8ml
- [CUY497P3](#) MS Stand Model Chamber Platinum Electrode
3mm gap L80xW3xH5mm 1.2ml
- [CUY497P4](#) MS Stand Model Chamber Platinum Electrode
4mm gap L80xW4xH5mm 1.6ml
- [CUY499L](#) Container for MS Stand Model Chamber,
Longer Type (Glass at bottom, For microscope
observation)

MS Stand Model Chamber SS Electrodes & Container



[CUY498-3](#) [CUY498-4](#) [CUY499](#)

- [CUY498-3](#) MS Stand Model Chamber SS Electrode 3mm
gap L50xW3xH10mm 1.5ml
- [CUY498-4](#) MS Stand Model Chamber SS Electrode 4mm
gap L50xW4xH10mm 2.0ml
- [CUY499](#) Container for MS Stand Model Chamber (Glass
at bottom, For microscope observation)

MS Gold Wire Electrode



[CUY500G0.5](#) [CUY500G1](#) [CUY500G2](#)

- [CUY500G0.5](#) MS Gold Wire Electrode 0.5mm gap
12.5µl 3/pkg.
- [CUY500G1](#) MS Gold Wire Electrode 1mm gap 25µl
3/pkg.
- MS Gold Wire Electrode 2mm gap 50µl

MS Gold Block Electrodes




[CUY501G1](#) [CUY501G2](#)

- [CUY501G1](#) MS Gold Block Electrodes 1mm gap 10µl
- [CUY501G2](#) MS Gold Block Electrodes 2mm gap 20µl

Electrodes for Nuclear Transfer (Animal Cloning)

Petridish Platinum Wire Electrodes for Activation




CUY5000P0.5

CUY5000P1

CUY5000P0.5 Petridish Platinum Wire Electrodes for Activation 0.5mm gap

CUY5000P1 Petridish Platinum Wire Electrodes for Activation 1mm gap

Electrodes for Nuclear Transfer Cloning (Chopstick Type) 10 μ m diameter, 25 μ m diameter




CUY5100-10

CUY5100-25

CUY5100-10 Electrode for Nuclear Transfer Cloning, 5mm-Tip, 10 μ m diameter (Chopstick Type) 2/pkg.

CUY5100-25 Electrode for Nuclear Transfer Cloning, 5mm-Tip, 25 μ m diameter (Chopstick Type) 2/pkg.

Electrodes for Nuclear Transfer Cloning (Chopstick Type) 50 μ m diameter, 100 μ m diameter



CUY5100-50

CUY5100-100

CUY5100-50 Electrode for Nuclear Transfer Cloning, 5mm-Tip, 50 μ m diameter (Chopstick Type) 2/pkg.

CUY5100-100 Electrode for Nuclear Transfer Cloning, 5mm-Tip, 100 μ m diameter (Chopstick Type) 2/pkg.