

# Electrodes

## CUY610P1.5-1, P2-1, P2.5-1, P4-1

Parallel Fixed Needle 45 Degree-Shape Exposed Platinum Electrode

0.5mm diameter

CUY610P1.5-1: L2-1, 1mm-Tip, 1.5mm gap

CUY610P2-1: L2-1, 1mm-Tip, 2mm gap

CUY610P2.5-1: L2-1, 1mm-Tip, 2.5mm gap

CUY610P4-1: L2-1, 1mm-Tip, 4mm gap, Standard Type (Nakamura Type)

