

## RFP antibody [6G6]


Product code: 6g6

*Note: NEW storage conditions and storage buffer, HIGHER purity.*

<b>Description</b>	Mouse monoclonal antibody [6G6] to Red Fluorescent Proteins
<b>Product type</b>	Primary antibody
<b>Host</b>	Mouse
<b>Isotype</b>	IgG2c
<b>Target/ Specificity</b>	Tested on mRFP, mCherry, mPlum, mOrange, mRFPruby, DsRed, mScarlet, mScarlet-I, OFF, tdTomato
<b>Immunogen</b>	UniProtKB: Q9U6Y8
<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	6G6
<b>RRID</b>	AB_2631395
<b>Concentration</b>	1 mg/mL
<b>Purity</b>	Protein A or G purified
<b>Form</b>	Buffered aqueous solution
<b>Storage buffer</b>	PBS, 50% glycerol, preservative: 0.09% sodium azide.
<b>Size</b>	20 µL; 150 µL
<b>Storage instructions</b>	Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.
<b>Applications</b>	Western blot: recommended starting dilution 1:2,000 The concentration of the antibody can vary. The optimal dilution should be determined by the user. A titration from 1:1,000 up to 1:5,000 is recommended. <i>Note: RFP antibody [6G6] does not work in immunofluorescence.</i>

### Tested applications

#### Western blotting



*Only for research applications, not for diagnostic or therapeutic use.*

*ChromoTek is a registered trademark of ChromoTek GmbH, part of Proteintech Group. Nanobody is a registered trademark of Ablynx, a Sanofi company. Alexa Fluor is a registered trademark of Life Technologies Corporation, a part of Thermo Fisher Scientific Inc. Other suppliers' products may be trademarks or registered trademarks of the corresponding supplier each. Statements on other suppliers' products are given according to our best knowledge.*